

# MR FIRST AID



24 - 8TH AVENUE EDENVALE, JOHANNESBURG, GAUTENG SOUTH AFRICA

TEL: 011-454 3570 FAX: 011-454 3472 WEB: [www.mrfirstaid.co.za](http://www.mrfirstaid.co.za) MAIL: [sales@mrfirstaid.co.za](mailto:sales@mrfirstaid.co.za)



## HEARTSAVE PAD

### **HeartSave SPEAKS YOUR LANGUAGE**

Clear voice prompts safely guide you through the rescue procedure and are available in a range of 40 different languages.

### **HeartSave IS FAST**

The device turns on automatically as soon as the cover has been removed and rhythm analysis starts immediately, after fixing the electrodes (SavePads).

### **HeartSave IS RELIABLE**

Clear, intuitive user interface, routine self-tests, protection against splash water, dust (IP55) and salt fog, extremely robust.

### **HeartSave IS SAFE**

**Thanks to sophisticated CCD technology from Primedic**

Defibrillation can be individually adjusted and is safe for the heart muscle.

### **Thanks to efficient data management**

All relevant data (ECG, buttons pressed, defibrillation's) is recorded and is immediately visible using the read-out software (AED, AS, optional with PAD).

### **Thanks to its update capability**

Software updates are easy to install via the SaveCard (memory card).

### **Thanks to clear status information**

The active status display provides up-to-date information on the usability of the device and battery.

### **Thanks to a long battery life**

High-performance lithium batteries last 3 (PAD / standby) to 6 (AED, AS / standby) years.

# HEARTSAVE PAD



SUCCESSFUL DEFIBRILLATION -

100-70%

70-50%

<50%

IF DEFIBRILLATION WITHIN  
3 MINUTES



IF DEFIBRILLATION WITHIN  
3-5 MINUTES



IF DEFIBRILLATION WITHIN  
5-10 MINUTES



## RESUSCITATION GUIDELINE ACCORDING TO ERC (EUROPEAN RESUSCITATION COUNCIL)

<p>1</p> <p>UNCONSCIOUS PERSON</p>	<p>2</p> <p>MAKE AN EMERGENCY CALL</p>	<p>3</p> <p>30 COMPRESSIONS : 2 RESCUE BREATHS</p>	<p>4</p> <p>FETCH / SEND FOR AED</p>	<p>5</p> <p>30 COMPRESSIONS : 2 RESCUE BREATHS</p> <p>CONTINUE 30:2 UNTIL AED IS AVAILABLE</p>
<p>6</p> <p>OPEN AED AND PLACE ELECTRODES DIRECTLY ON THE CHEST</p> <p>FOLLOW AED INSTRUCTIONS</p>	<p>7</p> <p>SHOCK ADVISED</p> <p>1 SHOCK</p> <p>2 MINUTES OF 30:2 RESUSCITATION STRAIGHT AWAY</p>	<p>7</p> <p>SHOCK NOT ADVISED</p> <p>2 MINUTES OF 30:2 RESUSCITATION STRAIGHT AWAY</p>	<p>8</p> <p>CONTINUE UNTIL EMERGENCY HELP ARRIVES!</p>	

